



AUTO TRANSFUSOR

Every drop counts...

Collect – anticoagulate – filter – online delivery
Collect – anticoagulate – filter – store – bedside delivery





AUTO TRANSFUSOR

AN ALL INTEGRAL BLOOD RECOVERY CUM
MANAGEMENT SYSTEM

PRODUCT FEATURES:

- Modularised design enables fast and safe handling of emergencies
- Offers affordable yet state of art solutions for harvesting blood components with down stream therapy specific filters
- Easy to setup and Easy to dismantle
- User-friendly and detachable filtration & storage modules
- Click-O-lock connector for blood bag attachment & detachment
- Contamination free handling of blood from surgical site to filtration, collection & storage
- Vacuum assisted transfer enables fast filtration & collection
- Facilitates on-line mixing of anticoagulant, if required
- Biocompatible Materials meets USP Class VI standards
- Crystal clear parts enable best viewing

TECHNICAL SPECIFICATIONS:

- Vacuum – 110 mmHg
- Micro-aggregate filter – 37 Micron
- Volume – 1100 ml



E-mail: enquiry@spictra.com
Web: www.spictra.com

Designed, Developed, Manufactured & Marketed by:

MEDIHUB SCIENCETEC PRIVATE LIMITED
Plot No:4, NH-7, MMDA Industrial Estate
Maraimalai Nagar - 603 209, Tamilnadu, India
Tel: +91 44 4740 0478, Fax: +91 44 4740 0484

In Collaboration with Direction Et Priorities, France

CUSTOMER CARE - **1800 102 2032**