



BED SIDE VACUUM CONTROLLER

*Touching Hearts.
Saving lives...*





BED SIDE VACUUM CONTROLLER


The Spictra bed side vacuum controller interfaces centralised vacuum line with bed side chest drainage system to provide safe control of the vacuum levels.

PRODUCT FEATURES

- Models: 1000ml and 750 ml capacity.
- Bio Compatible Materials.
- Rigid Body assures ZERO-COLLAPSABILITY during vacuum.
- Eco-Friendly, 100% disposable upon incineration.
- Compact and easy to install.

SPECIFICATIONS

- Housing : Styrene
- Vacuum Port : 1/4"
- CDS Port : 1/4"
- Atmosphere Line : 1/4"



Designed, Developed, Manufactured & Marketed by:

MEDIHUB SCIENCETEC PRIVATE LIMITED
Plot No:4, NH-7, MMDA Industrial Estate
Maraimalai Nagar - 603 209, Tamilnadu, India
Tel: +91 44 4740 0478, Fax: +91 44 4740 0484

E-mail: enquiry@spictra.com
Web: www.spictra.com

CUSTOMER CARE - **1800 102 2032**