

40µ BLOOD TRANSFUSION FILTER

Every drop counts...




ISO 9001:2015 Certified for Design, Manufacturing and Supply
ISO 13485:2016 Certified Medical Devices Platform
CE Marked Medical Devices



40 μ BLOOD TRANSFUSION FILTER

PRODUCT FEATURES:

- Depth filtration maximizes flow & throughput
- Traps at surface the Microaggregate, clots & particulate debris and also through the media
- Hard spike allows easy penetration of Blood Bags
- Meets emergency situations by handling multiple units of blood with high flow rates
- The integrated administration set quickens setup
- Very low held-up volume maximizes recovery



- 40 Micron Precision Filtering assures total safety
- Multi Stage Filtration design gives high throughput
- Held-up volume is just 4 ml



E-mail: enquiry@spictra.com
Web: www.spictra.com

Designed, Developed, Manufactured & Marketed by:

MEDIHUB SCIENCETEC PRIVATE LIMITED
Plot No:4, NH-7, MMDA Industrial Estate
Maraimalai Nagar - 603 209, Tamilnadu, India
Tel: +91 44 4740 0478, Fax: +91 44 4740 0484

CUSTOMER CARE - **1800 102 2032**