



ARTERIAL BLOOD FILTER AND BUBBLE TRAPS for CPB


Touching Hearts. Saving lives...



ISO 9001:2015 Certified for Design, Manufacturing and Supply
ISO 13485:2016 Certified Medical Devices Platform
CE Marked Medical Devices



ARTERIAL BLOOD FILTER AND BUBBLE TRAPS for CPB



PRODUCT FEATURES:

- Double Barrier assures total safety
- Long flow path maximizes reaction time
- Circular-O-Elliptical flow eliminates eddies
- Biocompatible Materials meets USP Class VI Standards
- Crystal Clear parts enables best viewing
-All these with lowest priming volume

Device	Model	Priming Volume, ML	Max blood flow rate, LPM	Filter size
Bubble Traps	Adult,	200	8	200 μ , Polyester
	Adult LP	120	8	200 μ , Polyester
	Infant & Paediatric	25	3	200 μ , Polyester
Arterial Blood Filter	Adult	120	8	200 μ , Polyester 60 μ , Polyolefin



Designed, Developed, Manufactured & Marketed by:

MEDIHUB SCIENCETEC PRIVATE LIMITED
Plot No:4, NH-7, MMDA Industrial Estate
Maraimalai Nagar - 603 209, Tamilnadu, India
Tel: +91 44 4740 0478, Fax: +91 44 4740 0484

E-mail: enquiry@spictra.com
Web: www.spictra.com

CUSTOMER CARE - **1800 102 2032**