



# CONCENTRIC NEEDLE ELECTRODE FOR EMG

*Precise Technology...  
Precise Diagnosis*





# CONCENTRIC NEEDLE ELECTRODE FOR EMG

## PRODUCT FEATURES:

- Electrode variants – 0.5mm diameter – 25mm, 40mm & 50mm lengths
- Meets clinical requirements for different muscles and patient profile ranging from paediatric to adults
- 0.07 mm<sup>2</sup> pickup area and cable assures noise free recordings with matching impedance
- 3 pin or 5 pin cables give connectivity with all machines
- Spictra Concentric Needle Electrodes are sterilized by EtO
- All needles are meant for single use only
- Destroy after use



E-mail: [enquiry@spictra.com](mailto:enquiry@spictra.com)  
Web: [www.spictra.com](http://www.spictra.com)

Designed, Developed, Manufactured & Marketed by:

**MEDIHUB SCIENCETEC PRIVATE LIMITED**  
Plot No:4, NH-7, MMDA Industrial Estate  
Maraimalai Nagar - 603 209, Tamilnadu, India  
Tel: +91 44 4740 0478, Fax: +91 44 4740 0484

In Collobration with SCTIMST, Trivandrum

CUSTOMER CARE - **1800 102 2032**