



HEMOCONCENTRATOR




Touching Hearts. Saving lives...



ISO 9001:2015 Certified for Design, Manufacturing and Supply
ISO 13485:2016 Certified Medical Devices Platform
CE Marked Medical Devices



HEMOCONCENTRATOR



Hematocrit: $32 \pm 2\%$
 Protein: $60 \pm 5\text{g/L}$
 Temperature: $37 \pm 1^\circ\text{C}$
 Bicarbonate: $20 \pm 2 \text{ mg/dL}$

Fibre area	Priming Volume	Max Flow rate	Max TMP
0.7m2	46ml	500 ml/min	500 mm Hg
0.3m2	28ml	500 ml/min	500 mm Hg



E-mail: enquiry@spictra.com
 Web: www.spictra.com

Designed, Developed, Manufactured & Marketed by:

MEDIHUB SCIENCETEC PRIVATE LIMITED
 Plot No:4, NH-7, MMDA Industrial Estate
 Maraimalai Nagar - 603 209, Tamilnadu, India
 Tel: +91 44 4740 0478, Fax: +91 44 4740 0484

In Collobration with SCTIMST, Trivandrum

CUSTOMER CARE - **1800 102 2032**