



MYOCARDIAL PROTECTION SYSTEM

*Touching Hearts.
Saving lives...*





ISO 9001:2015 Certified for Design, Manufacturing and Supply
ISO 13485:2016 Certified Medical Devices Platform
CE Marked Medical Devices



MYOCARDIAL PROTECTION SYSTEM

PRODUCT FEATURES

- Adult and Paediatric delivery system configurations
- Flow engineering with Axial, non turbulent, eddy free design, prevents bubble formation and eliminates need for downstream filters
- Shortest blood flow path over smooth surfaces .. yet best cooling efficiency
- Core flow without thinning of blood, reduces stress
- No held up volume due to top-down reversible design feature
- Mounting direction aspecific
- Priming volume is just 40 ml
- Mixing proportions of blood & Cardioplegia 1:1, 2:1, 4:1



Designed, Developed, Manufactured & Marketed by:

MEDIHUB SCIENCETEC PRIVATE LIMITED
Plot No:4, NH-7, MMDA Industrial Estate
Maraimalai Nagar - 603 209, Tamilnadu, India
Tel: +91 44 4740 0478, Fax: +91 44 4740 0484

E-mail: enquiry@spictra.com
Web: www.spictra.com

CUSTOMER CARE - **1800 102 2032**