

CARDIOPLEGIA DELIVERY SYSTEM

for Crystalloid Cardioplegia &
Blood Cardioplegia

*Touching Hearts.
Saving lives...*




ISO 9001:2015 Certified for Design, Manufacturing and Supply
ISO 13485:2016 Certified Medical Devices Platform
CE Marked Medical Devices





CARDIOPLEGIA DELIVERY SYSTEM

UNIQUE "HELICAL – COUNTER CURRENT – DOWN FLOW" ACHIEVES...

- Gentle mixing of crystalloid & blood
- Lowest contact area .. smooth surface .. yet fastest cooling
- Elimination of air bubbles
- Integrated course filter



Helical flow enhances path length and hence heat transfer



E-mail: enquiry@spictra.com
Web: www.spictra.com

Designed, Developed, Manufactured & Marketed by:

MEDIHUB SCIENCETEC PRIVATE LIMITED
Plot No:4, NH-7, MMDA Industrial Estate
Maraimalai Nagar - 603 209, Tamilnadu, India
Tel: +91 44 4740 0478, Fax: +91 44 4740 0484

CUSTOMER CARE - **1800 102 2032**