



CHEST / THORACIC DRAINAGE SYSTEM

PRODUCT FEATURES

- Closed System with separate filling / vacuum port
- Ensures that sterility is not lost while filling water in the product
- Pre-connected Patient Line
- Ensures contamination-free setup
- Over-flow design
- Offers précised measurements in calibrated graduations
- Absolutely Transparent
- Transparent crystal clear, round reservoir offers 360° viewing of contents
- Rigid Reservoir assures ZERO-COLLAPSABILITY during vacuum
- Sloped bottom, yet standing flat
- Ensures safe water seal

ADVANTAGES

- Over flow model: Offers higher height for even low volumes enhancing drain diagnosis during initial post-operative period
- One device for both Adult & Pediatric
- Total protection against ascending infection
- Multiple Models – Multiple advantages



E-mail: enquiry@spictra.com
Web: www.spictra.com

Designed, Developed, Manufactured & Marketed by:

MEDIHUB SCIENCETEC PRIVATE LIMITED
Plot No:4, NH-7, MMDA Industrial Estate
Maraimalai Nagar - 603 209, Tamilnadu, India
Tel: +91 44 4740 0478, Fax: +91 44 4740 0484

CUSTOMER CARE - **1800 102 2032**